



# Thermalox TOC-TN

The modern alternative to COD and TKN

# THE SERCON THERMALOX

## Measuring Total Organic Carbon

TOC (Total Organic Carbon) is a measurement widely used in a number of industries to assess the organic contamination of water. Applications include Environmental Research, Industrial Process Monitoring, and Quality Assurance.

In wastewater TOC concentrations of many thousands of ppm are often found. At the other extreme, low ppb TOC levels must be measured in the Ultra-Pure water used by Pharmaceutical, Semiconductor and the Power Generation industries.

There are several Oxidation methods which have been historically utilised to measure TOC. The Sercon Thermalox uses the Thermal Catalytic Oxidation technique which gives excellent performance for a wide range of sample types and applications.

| OXIDATION METHOD            | SPECIFICATIONS  | DISADVANTAGES  |
|-----------------------------|---|--|
| Wet Chemical                | Inexpensive and robust.   | Inefficient oxidation. Poor with particulate or hard to oxidise compounds. Slow speed of analysis. |
| UV Persulfate               | Good Precision, Low detection limits. Ideal for DOC in pure waters. | Poor performance with particulate or humic material (i.e. natural waters)                          |
| Thermal Catalytic Oxidation | Near perfect recovery (oxidation). Fast.                            | Higher ownership costs. Trace catalyst contamination means sensitivity is compromised              |



# THE SERCON THERMALOX TOC/TN<sub>B</sub>

---

Complete Oxidation and detection are critical for measuring TOC and TN<sub>B</sub>. Working with one of the world's leading NDIR analyser manufacturers, has dramatically increased the sensitivity of the Sercon Thermalox CO<sub>2</sub> detector. The Sercon Thermalox detector is now over twenty times more sensitive than its leading rival.

---

Additionally, the combination of catalytic combustion and unique furnace design allows efficient analysis of organic compounds whilst virtually eliminate trace carbon contamination within the reactor and catalyst.

High Sensitivity allows lower injection volumes to be used leading to less catalyst contamination, more accuracy, and faster measurements. It also provides a dramatic improvement in the detection limit for this method.

Other methods employed to reach these low detection limits result in reduced performance in terms of increased maintenance and longer analysis times.

The Sercon Thermalox NO<sub>x</sub> detector employs a number of features to ensure excellent sensitivity and stability. An innovative NO<sub>x</sub> → NO reduction furnace is used which remains efficient over an extended life. The detection reaction for TN is heated under vacuum using precise temperature and vacuum controls guaranteeing stability and optimised sensitivity.

One of the benefits of the unique Sercon Thermalox detectors is the long-term stability of the calibration, giving the user more time on analysis and less time on set up.

These improvements combined with, the sophisticated but user-friendly software and excellent build quality, highlight the Sercon Thermalox as the most cost-effective, powerful and robust tool for the analysis of TOC and TN<sub>B</sub> anywhere in the world.

## System Highlights:

- Catalytic thermal oxidation. The only way to measure complex samples containing particulate or difficult to oxidise materials
- Standard Deviation better than 20ppb on low level Carbon
- Standard Deviation better than 10ppb on low level Nitrogen
- Upper Range Limit greater than 50,000ppm for TOC
- Detection Limit of better than 40ppb for TOC and 10ppb for TN<sub>B</sub>
- Short Analysis time of less than two minutes
- Complete recovery, including Suspended Solid fraction
- Handles salts and particulate easily without sacrificing detection limits or requiring specific combustion tubes
- Capable of analysing small samples as standard.
- Holds up to 88 samples
- Complete washing between samples ensures no carry over
- Totally software driven from a Windows™ based platform
- Automatic preparation of calibration standards
- Versatile to your analysis – you can have standalone TOC, standalone TN<sub>B</sub> or combination TOC/TN<sub>B</sub>
- The Sercon Thermalox is so well automated that you can mix high and low level samples on the same vial rack

# THE SERCON THERMALOX

## Features included as standard

The Sercon Thermalox instrument uses an optimised XYZ autosampler. This gives maximum flexibility and allows a high degree of automation, whilst maintaining robust performance.


The standard vial rack holds up to 88 glass or disposable plastic vials, with or without lids. Sample racks can be chilled, and optional racks for special requirements are available.

Calibrants can be made up automatically; sparging, dilutions and sample agitation can easily be carried out all under software control.

By using a special aspiration step, samples containing particulate can be automatically agitated, ensuring that a representative sample aliquot is aspirated from the vial.

- Both TOC by subtraction and by pre-stripping of inorganic carbon methods are fitted as standard
- Carrier gas generation fitted as standard on TOC – no need for bottled gases
- Automatically prepares its own calibrants – reducing possibility of error and demands on user time
- Peltier cooler for condensate removal - no downtime or costs involved with changing desiccants
- Precise electronic mass flow control of the carrier gas – stable calibrations and results
- Highly sensitive detectors – prolongs catalyst life, lowers detection limit
- Stand-by Mode - limits carrier gas consumption and prolongs oven life improving cost of ownership
- Specially manufactured Carbon-free catalyst - virtually no blank carbon peak
- Direct Injection method - eliminates blockage or carryover contamination
- **Sample sparging and POC analysis method included as standard**

## Windows Platform Software



- Intuitive functionality
- User friendly GUI
- Secure to FDA21 CFR11
- Robust
- Tailored to suit you



# THE SERCON THERMALOX

## Total Nitrogen


Total Nitrogen (TN) determination may be added to The Sercon Thermalox TOC analyser to give TOC and TN measurements from the same vial. Sercon also offer The Sercon Thermalox TNb analyser where only this analyte is required.

In the same way as TC measurements are performed, a sample aliquot is injected into a catalytic furnace. Carrier gas (oxygen) sweeps the oxidant, in this case NOx gases, through to the NOx detector.

The detector utilizes the chemiluminescent reaction between NO and O<sub>3</sub> to determine the NO concentration. As well as a reduction step to account for NO<sub>2</sub> species present. The importance of this reduction step is illustrated by the different NO/NO<sub>2</sub> splits of the various compounds shown below.

### Oxidation Ratios NO/NO<sub>2</sub>

■ NO ■ NO<sub>2</sub>




The unique features of the chemiluminescence detector gives The Sercon Thermalox TN analyser its remarkable sensitivity, at least an order of magnitude more than other available systems, by calibration, TNb is determined from the amount of NO measured.

As can be seen from the graph below, temperature, concentration and injection volume have a critical effect on the conversion of R-N to NO. For most applications, the TC catalytic furnace can be modified to allow it to perform both the TC and TN measurements. However, for some higher level TN applications – up to 200mg/l – or highly particulated samples, a specialised TN catalyst is added.

This overcomes the quenching effects experienced on many TN applications where larger sample volumes are needed or higher concentrations of TN are measured - for example effluent samples. Large bore needles may be used to aspirate particulated samples representatively.

### 80 µg/l Injection, Creatinine



# THE SERCON THERMALOX TOC/TN<sub>B</sub>

## RANGEABILITY

Calibration – from two to 7 points may be automatically defined in the method set-up. The software then uses regression mathematics to generate a best fit curve. This may be a linear or an n-order polynomial function to give unrivalled rangeability.

## AUTOMATION

Users can select a method which includes automatically preparing calibrants from a stock solution, calibration, acid stripping and measurement of TOC and TN all in one operation. The principals of FDA21 CFR11 are applied to ensure that the data is tamperproof and auditable.


## SENSITIVITY AND PRECISION

It's very hard to beat the sensitivity and precision of The Sercon Thermalox range of elemental analysers. The following results were obtained on The Sercon Thermalox TOC/TN\*.

\*The samples were prepared by adding each of the chemicals in the nominal concentrations to deionised water. The digestion efficiency represents the ratio of nitrogen measured compared to the amount of chemically bound nitrogen contained in the sample. Oven Temperature was 950°C. Lower oven temperatures give a wider spread of recoveries at higher concentrations but have some advantages with some matrices at low levels.

| TOC/TN <sub>0</sub> TEST | TC µg/l KHP | TC mg/l KHP | TC µg/l KHP | TN µg/l NH <sub>4</sub> Cl | TN mg/l NH <sub>4</sub> Cl | TN mg/l NH <sub>4</sub> Cl |
|--------------------------|-------------|-------------|-------------|----------------------------|----------------------------|----------------------------|
| Concentration            | 145         | 1.97        | 10,108      | 101                        | 1.03                       | 201.6                      |
|                          | 138         | 1.99        | 10,023      | 98                         | 1.02                       | 200.8                      |
|                          | 143         | 2.00        | 10,098      | 98                         | 0.98                       | 199.4                      |
| Mean                     | 142         | 1.99        | 10,076      | 99                         | 1.01                       | 200.6                      |
| Nominal Concentration    | 140         | 2.00        | 10,000      | 100                        | 1.00                       | 200.00                     |
| Concentration            |             |             |             |                            |                            |                            |
| SD                       | 2.9µg/l     | 0.01mg/l    | 37.93mg/l   | 1.25µg/l                   | 0.02mg/l                   | 0.91mg/l                   |
| CV                       | 2.0%        | 0.5%        | 0.38%       | 1.25%                      | 2.0%                       | 0.45%                      |

| TN Digestion Efficiency 40µl Injection | Creatinine | NH <sub>4</sub> Cl (Calibrant) | NaNO <sub>3</sub> |
|--|------------|--------------------------------|-------------------|
| Concentration (mg/lN)                  | 98.8       | 100.3                          | 102.0             |
|  | 98.2       | 99.8                           | 101.8             |
|  | 98.7       | 100.1                          | 101.7             |
| Mean                                   | 98.6       | 100.10                         | 101.87            |
| Nominal Concentration                  | 100.00     | 100.00                         | 100.00            |
| SD (mg/l N)                            | 0.26       | 0.20                           | 0.13              |
| Digestion Efficiency                   | 98.6%      | 100.1%                         | 101.9%            |



# THE SERCON THERMALOX

|                              |   |
|------------------------------|---|
| Analyte                      | TC, TIC, TOC, NPOC, POC with options for measuring TN and TP  |
| Applications                 | Drinking water, pharmaceutical cleaning in place water, ground water, surface water, saline waters, domestic, and industrial wastewater   |
| Method                       | Total substance: Thermal catalytic oxidation at up to 1000°C. Inorganic substance: acid stripping or injection into low temperature acid reactor (120°C)  |
| Range                        | Carbon: Various from less than a range of 40-150µg/l to greater than a range of 100-40,000mg/l: Nitrogen: From 20-100 g/l to 1.0-200mg/l  |
| Detection limit              | 40 ppb for TOC, 10ppb of TN   |
| Auto Dilution                | 500:1   |
| Cycle time                   | Typically 120 – 180 seconds   |
| Precision                    | Standard deviation: ≤ 5% of full scale for ranges to 3mg/l; ≤ 3% of full scale for ranges up to 500mg/l and 2% of full scale for ranges up to 40,000mg/l  |
| Sample matrix                | Aqueous samples including those containing suspended solids and salts. Acids, slurries, sludges and organics by special application. A solids option is available   |
| Injection volume             | 3 to 250µl. Max 1,000µl ( option)<br>Needle diameters up to 1.0mm are available for particulate samples   |
| Sample injection             | Automatically from an 88 position vial tray, or manually  |
| Carrier gas                  | A carrier gas generator is fitted as standard for TOC. For TN, oxygen (grade 5.5 or better) is required; flow rate is 15 litres per hour or less  |
| Software                     | Windows™ base PC control English, multi language (option)   |
| Power supply                 | 230V AC ± 10%, 50Hz ± 1%; or 115V AC ± 10%, 60 Hz ± 1%  |
| Ambient temperature          | 5 – 35°C  |
| Dimension and weight         | 645W x 505D x 650H mm footprint (WDH) ; weight 50kg; (figures include sampler)  |
| Protection Class             | EMI class Euro 50081/50082  |
| TOC and TN method compliance | BS and DIN EN 1484 of 1997; ISO 8245, AOAC973.47, Standard Methods 5310B, Environmental Protection Agency (EPA) Methods 415.1, EPA SW-846 Method 9060A, ASTM D2579-85. TN according to DIN38409 part 27, ENV 12260 and ISO/TR11905-2. |

Sercon Limited, Unit 3B, Crewe Trade Park,  
Gateway, Crewe, Cheshire CW1 6JT UK  
Tel: +44 (0)1270 580008 Email: sales@sercongroup.com

[www.sercongroup.com](http://www.sercongroup.com)

 SerconLimited  @SerconLimited



ISO 9001:2015 Certified  
ISO 13485:2016 Certified

December 2020



COLLABORATE

CREATE

INNOVATE

Unit 3B Crewe TradePark,  
Crewegates,  
Crewe,  
Cheshire.  
CW1 6YY

---

t: +44 (0)1270 580008  
f: +44 (0)1270 252310  
e: [info@sercongroup.com](mailto:info@sercongroup.com)

---

[www.sercongroup.com](http://www.sercongroup.com)