



ONE AUTOSAMPLER FOR BOTH OILS AND AQUEOUS

switch between
oils and **coolants** testing
on demand

- Mix samples immediately prior to analysis for accurate results
- Unique dual-matrix design for fast switching between sample types
- Fast, smooth XYZ movement
- Accurate positioning to ± 0.4 mm
- No probe/stirrer alignment drift



Oils 7400 & Oils 7600 HOMOGENIZING AUTOSAMPLERS




Why sample mixing?

Used oil samples typically carry particulate material that can settle prior to analysis; this can lead to non-representative data. The Oils 7400 and Oils 7600 automation have been designed to resolve this problem.

Prior to analytical measurement, a stirring paddle mounted next to the sample probe efficiently mixes each sample and, with the sample probe, is cleaned at the rinse station. Sample mixing is configurable, via the dashboard software, to meet the challenge of more viscous oil samples.


Accurate and fast

Each sample is automatically homogenized prior to analysis. While the current sample is drawn, the next sample is being stirred to increase efficiency and save time.




Maintain sample integrity

An integrated drip cup travels with the sample probe and stirring paddle to capture any sample droplets, eliminating cross contamination while moving over the sample tubes.




Quickly switch samples

Two discrete liquid flow paths and a dual port rinse station allow the user to switch between oils and coolants testing on demand in as little as 5 minutes.





Filtered probe

The included stainless steel filtered probe helps filter out large particles, preserving the lifetime of downstream components. The material is robust and easy to clean.



Technical Specifications

Dimensions

		
	Oils 7400	Oils 7600
Height	49 cm (19.3 in)	49 cm (19.3 in)
Width	57 cm (22.4 in)	81 cm (31.9 in)
Depth	57.5 cm (22.7 in)	57.5 cm (22.7 in)
Weight	23 kg (50 lbs)	35 kg (77 lbs)

Capacities

Oils 7400: Up to 4 racks, up to 384 samples

Oils 7600: Up to 6 racks, up to 576 samples

Rack options

Bel-Art, CETAC/Bohdan, Janus

Hardware interfaces

RS-232 and USB

Power requirements

100-240 VAC, 47-63 Hz, 1.9 A

Minimum computer requirements

Microsoft Windows® 10 or later operating system

1 GB RAM

2 COM or USB port

Internet Explorer 6 or higher must be installed

for system to function properly

Warranty

2 year limited

Teledyne CETAC Technologies, 14306 Industrial Road, Omaha, NE 68144 USA
+1 402.733.2829 | +1 800.369.2822 | cetacsales@teledyne.com | www.teledynecetac.com

©2022 TELEDYNE TECHNOLOGIES Specifications are subject to change without notice. Other products and company names mentioned herein may be trademarks and/or registered trademarks.

Revised 2022 November 29