

Titramax 1000

Technical datasheet



Dimensions (W × H × D)	320 × 128 × 375 mm
Weight	8 kg
Maximum permissible load	5 kg
Drive	Shaded pole motor
Speed range	150 – 1,350 rpm
Acoustic pressure	< 50 dB (A)
Rated voltage	n.a.
Overvoltage category	II
Degree of pollution	2
Power input	25 W
Operating temperature	5 °C – 31 °C up to 80 % rel. humidity 32 °C – 40 °C up to 50 % rel. humidity (decreasing linearly)
Installation altitude	up to 2,000 m asl
Overheat protection	self-resetting
Speed control	analogue
Usable area	290 × 258 mm
Movement	circular vibrating
Orbit	1.5 mm
Fuse	M 1,25 A (230 V) or T 2,0 A (115 V)
Permissible mains voltage fluctuations	±10 %
Bottle dimensions	n.a.

Interfaces

n.a.

Protection Class