

Hei-Mix S

Technical datasheet



Dimensions (W × H × D)	126 × 80 × 140 mm
Usable surface heating plate	104 mm stirrer plate
Weight	1.1 kg
Maximum permissible load	6 kg
Drive	shaded pole motor
Speed range	0 – 2.200 rpm
Heating power	n.a.
Heating temperature range	n.a.
Heating control	n.a.
Temperature measurement accuracy	n.a.
Temperature measurement resolution	n.a.
Display	n.a.
Protection class (EN 60529)	IP 21
Acoustic pressure	< 70 db(A) (in accordance with IEC 61010)
Rated voltage	n.a.
Connection	n.a.
Protection Class	2
Overvoltage category	n.a.
Degree of pollution	n.a.
Power input	n.a.
EMC class	n.a.

Operating temperature	n.a.
Installation altitude	n.a.
Indication/measurement range	n.a.