

Hei-VAC Valve Industrial

Technical datasheet



Dimensions (W × H × D)	326 x 243 x 402 mm
Weight	17.3 kg
Protection class (EN 60529)	IP 40
Connection	Coolant: hose nozzle DN 6-8 mm
Operating temperature	+10 - +40°C
Motor protection	thermal cutout, manual reset
Storage temperature range	-10 - +60°C
Power	Electrical: 0.25 kW
Maximum pumping speed	3.4 / 3.8 m ³ /h
Ultimate vacuum (absolute)	without gas ballast: 1.5 mbar with gas ballast: 3 mbar
Maximum permissible pressure	1.1 bar
Permissible relative atmospheric moisture during operation (no condensation)	30% - 85%
Maximum no-load speed	1500 / 1800 min ⁻¹
Inlet	hose nozzle DN 8 mm
Outlet	hose nozzle DN 8-10 mm
Maximum permissible range of supply voltage (±10%)	100-115 V~ 50/60 Hz, 120 V~ 60 Hz / 200-230 V~ 50/60 Hz