

JLT Series

Flocculators

Multi-position digital Jar Testers with timer and illuminated back panel.

velp®

Driven by curiosity



Easy to use and Reliable

Digital Flocculators in 4 and 6 positions with timer and illuminated back panel to allow standard conditions to be met during Jar Testing in water treatment plants.

- Digitally selectable speed from 10 to 300 rpm
- Microprocessor controlled operating timer for unattended operation
- Stirring rods adjustable in height without tools
- Intuitive and bright digital display to precisely set stirring speed and timer
- Identical speed and time for all positions for comparable results
- Strong resistance to chemical and mechanical corrosion



3 YEARS WARRANTY

Instrument - Codes

JLT 4	100-240 V / 50-60 Hz	F105A0108
JLT 6	100-240 V / 50-60 Hz	F105A0109

Technical Data

	JLT 4	JLT 6
Construction material	Epoxy painted metal structure	Epoxy painted metal structure
Number of samples	4 positions	6 positions
Back panel	Disconnectable illuminated	Disconnectable illuminated
Stainless steel stirring rods	Adjustable in height by a selfblocking chuck	Adjustable in height by a selfblocking chuck
Monoselector	Same speed for each rod	Same speed for each rod
Electronic speed control	From 10 to 300 rpm	From 10 to 300 rpm
Speed setting interval	1 rpm	1 rpm
Microprocessor controlled timer	0÷999 min or 0 ÷99 hours (or continuous)	0÷999 min or 0 ÷99 hours (or continuous)
Power	11 W	19 W
Dimensions (WxHxD)	645x347x260 mm (25.4x13.7x10.2 in)	935x347x260 mm (36.8x13.7x10.2 in)
Weight	13 kg (28.6 lb)	17 kg (37.4 lb)



Headquarters
Velp Scientifica Srl
Via Stazione 16
20865 Usmate (MB)
Italia
T +39 039 628811
velpitalia@velp.com

LATAM
velplatam@velp.com

USA
Velp Scientifica Inc
40, Burt Drive, Unit #1,
Deer Park
NY 11729 - U.S.
T +1 631 573 6002
velpusa@velp.com

SEA & Pacific
velpsea@velp.com

China
Velp China Co. Ltd.
Room 828, Building 1,
778 Jinji Road, Pudong,
Shanghai, China
Tel. +8621 34500630
velpchina@velp.com

India
Velp India Pvt. Ltd.
Office 323-324, Building 2
Raheja Tesla, MIDC, Turbhe,
400705, Navi Mumbai, India
velpindia@velp.com

Velp Official Partner

Rev 3.02.2026

