




SP 311

Peristaltic Pumps

Peristaltic pumps ideal for transferring gas, liquids, mixtures and emulsions in various industries.



SP 311 Peristaltic Pumps

Easy to use and Reliable

- Minimum maintenance (just replace the tube)
- Easy cleaning – no contact of liquid and mechanical parts
- Tubing (2 diameters available) resists most chemicals
- Electronic speed control
- Pressure: 10 metres water column

Recommended for:

- Transfer of liquids, gases, mixtures, emulsions, etc.
- Taking samples and continuous dosing
- Analysis automation



INSTRUMENT - CODE

Instrument	Flows (ml/min.)	Tube dimensions (mm)	Code
SP 311/2	2 - 35	4x7 Tube dimensions (mm)	F40200002
SP 311/6	6 - 135	4x7 Tube dimensions (mm)	F40200006
SP 311/12	20 - 240	4x7 Tube dimensions (mm)	F40200012
SP 311/60	125 - 1000	5x8 Tube dimensions (mm)	F40200060

TECHNICAL DATA

Construction material	Epoxy painted stainless steel structure
Motor	24 V CC with gear reduction
Power	20 W
Power supply	230 V / 50 - 60 Hz
Weight	1,8 Kg (3.9 lb)
Dimensions (WxHxD)	145x95x190 mm (5.7x3.7x7.5 in)



We reserve the right to make technical alterations
We do not assume liability for errors in printing, typing or transmission

VELP Official Partner

DESIGNED AND MANUFACTURED IN ITALY



ITALY - HQ
Via Stazione 16
20865 Usmate (MB) Italy
Tel. +39 039 628811
velpitalia@velp.com

INDIA
velpindia@velp.com

USA
40, Burt Drive, Unit #1, Deer Park
NY 11729 - U.S.
Tel. +1 631 573 6002
velpusa@velp.com

LATAM
velplatam@velp.com

CHINA
Room 828, Building 1, No.778
Jinji Road, Pudong New Area,
Shanghai, China
Tel. +8621 34500630
velpchina@velp.com