

CLASSIC

Vortex Mixer

Compact vortex mixer suitable for different mixing requirements, with touch or continuous mode.

velp®

Driven by curiosity



Compact, Robust and Convenient

Compact vortex mixer suitable for different mixing requirements, with touch or continuous mode. It is highly resistant to chemicals and has 3 anti-slip rubber feet for stability.

- Automatic start in touch mode or continuous operation
- Variable speed operation, electronically controlled up to 3000 rpm
- Optimum resistance to chemicals from technopolymer casing
- Excellent stability with zinc alloy base
- IP 42 protection against particles and liquids



3 YEARS WARRANTY

Instrument - Codes

Classic

100-240 V / 50-60 Hz

F202A0173

Technical Data

Construction material	Zinc alloy and technopolymer
Support system	3 anti-sliding feet
Operational mode	Touch / Continuous
Speed setting	Analog
Speed control	Electronic
Stirring speed	0 - 3000 rpm
Orbital diameter	4,5 mm
Type of movement	Orbital
Electronic protection degree CEI EN 60529	IP 42
Power	15 W
Weight	2.2 kg (4.9 lb)
Dimensions (WxHxD)	180x70x220 mm (7.1x2.8x8.7 in)



Headquarters
Velp Scientifica Srl
Via Stazione 16
20865 Usmate (MB)
Italia
T +39 039 628811
velpitalia@velp.com

LATAM
velplatam@velp.com

USA
Velp Scientific Inc
40, Burt Drive, Unit #1,
Deer Park
NY 11729 - U.S.
T +1 631 573 6002
velpusa@velp.com

SEA & Pacific
velpsea@velp.com

China
Velp China Co. Ltd.
Room 828, Building 1,
778 Jinji Road, Pudong,
Shanghai, China
Tel. +8621 34500630
velpchina@velp.com

India
Velp India Pvt. Ltd.
Office 323-324, Building 2
Raheja Tesla, MIDC, Turbhe,
400705, Navi Mumbai, India
velpindia@velp.com

Velp Official Partner

Rev 3.02.2026

