

RC 5

Heating Plates

Single-position heating plate with aluminum top.

velp[®]

Driven by curiosity



Excellent Heating Performance

The RC 5 is equipped with a 135 mm diameter aluminum alloy top plate ensuring a great temperature homogeneity and optimum heat transfer across the entire surface.

Designed to guarantee ease of use, the longest lifespan and excellent resistance to chemicals.

- Maximum temperature 310 °C
- Accepts flasks of up to 15 L volume
- Hot plate warning indicators
- IP 42 certified



3 YEARS WARRANTY

Instrument - Codes

RC 5	230 V / 50-60 Hz	F20700610
RC 5	115 V / 50-60 Hz	F20700610

Technical Data

Housing material	Epoxy painted aluminum structure
Plate material	Aluminum alloy
Plate dimension	135 mm
Temperature range	Room temp. - 310 °C
Hot top warning system	illuminated icon (when the temperature is above 50 °c)
Power input	630 w
Permissible ambient temperature	5 - 40 °C
Protection class	IP 42
Weight	1.4 kg (3.1 lb)
Dimensions (WxHxD)	160x85x270 mm (6,3x3,3x10,6 in)



Headquarters
Velp Scientifica Srl
Via Stazione 16
20865 Usmate (MB)
Italia
T +39 039 628811
velpitalia@velp.com

LATAM
velplatam@velp.com

USA
Velp Scientific Inc
40, Burt Drive, Unit #1,
Deer Park
NY 11729 - U.S.
T +1 631 573 6002
velpusa@velp.com

SEA & Pacific
velpsea@velp.com

China
Velp China Co. Ltd.
Room 828, Building 1,
778 Jinji Road, Pudong,
Shanghai, China
Tel. +8621 34500630
velpchina@velp.com

India
Velp India Pvt. Ltd.
Office 323-324, Building 2
Raheja Tesla, MIDC, Turbhe,
400705, Navi Mumbai, India
velpindia@velp.com

Velp Official Partner

Rev 3.02.2026

