

REC 10 Digital Heating Plates

Digital heating plate with ceramic with a large surface
for demanding heating tasks.

velp®

Driven by curiosity



Large Surface Heating Plate

The REC 10 Digital provides a 10-inch spacious plate for demanding heating tasks, ensuring great temperature homogeneity and optimum heat transfer.

- Designed to guarantee ease of use and excellent resistance to chemicals
- Digital display for precise temperature settings up to 550 °C
- A "Hot Surface" warning message is displayed when the temperature exceeds 50 °C, even when inactive
- Accepts flasks of up to 25 L volume
- Precise thermoregulation of the medium with Pt100/Pt1000 probe (± 1 °C) and VTF Digital Thermoregulator (± 0.5 °C)
- Designed for the longest lifespan and guaranteed by 3 years warranty with online product registration



3 YEARS WARRANTY

Instrument - Codes

REC 10 Digital	230 V / 50-60 Hz	F20700693
REC 10 Digital	115 V / 60 Hz	F20710693
REC 10 Digital	100 V / 50-60 Hz	F20780693

Technical Data

User interface	Digital
Housing material	Aluminum body - technopolymer enclosure
Plate material	Ceramic
Plate dimension	260x260 mm; 10x10 in
Temperature range	Room temp - 550 °c
Heating control	Digital
External temperature sensor connection	Pt100 Pt1000 VTF
Thermoregulation range	Room temp. - 350 °c
Thermoregulation accuracy	± 1 °C (Pt100; Pt1000) ± 0.5 °C (VTF)
Hot plate warning system	On display (when the temperature is above 50 °c)
Power input	1510 W (230 v / 50-60 hz) 1500 W (115 v / 60 hz) 1300 W (100 v / 50-60 hz)
Weight	5.3 kg 11.6 lb
Dimensions (WxHxD)	280x100x420 mm 11.02x3.93x16.53 in



Headquarters
Velp Scientifica Srl
Via Stazione 16
20865 Usmate (MB)
Italia
T +39 039 628811
velpitalia@velp.com

LATAM
velplatam@velp.com

USA
Velp Scientific Inc
40, Burt Drive, Unit #1,
Deer Park
NY 11729 - U.S.
T +1 631 573 6002
velpusa@velp.com

SEA & Pacific
velpsea@velp.com

China
Velp China Co. Ltd.
Room 828, Building 1,
778 Jinji Road, Pudong,
Shanghai, China
Tel. +8621 34500630
velpchina@velp.com

India
Velp India Pvt. Ltd.
Office 323-324, Building 2
Raheja Tesla, MIDC, Turbhe,
400705, Navi Mumbai, India
velpindia@velp.com

Velp Official Partner

Rev 3.02.2026

